





MILESTONE
H E L P I N G
P A T I E N T S



PRESTO PRO


**Processor/Stainer for
High Quality Frozen Sections**




Automated Stainer

PRESTO PRO

PrestoPRO is an automated benchtop stainer designed to process and stain frozen section slides, providing exceptional H&E staining quality and reproducibility. By automating the staining process, PrestoPRO improves lab efficiency by freeing up the technician's time while simultaneously producing consistently high-quality results for the patient. The PrestoPRO system is a great addition to any frozen section lab, as all staining protocols are customizable and can be optimized to fit the customer's specifications and needs.




INCOMPARABLE RESULTS



Breast tissue processed with
Milestone's PrestoPRO+FineFIX Solution



Combine PrestoPRO with Milestone's PrestoCHILL
for superior quality frozen section slides



AUTOMATION:

Automate time in each station, eliminating staining variability



STANDARDIZATION:

Consistent staining results, no matter the operator



EFFICIENCY:

Improve workflow TAT in continuous batch mode



SAFETY:

Fully vented and enclosed chamber for fume control



FLEXIBILITY:

Set and modify protocols according to your needs

KEY BENEFITS

HIGH QUALITY SLIDES

- Frozen section slide quality that rivals FFPE slides
- Improved morphology
- Fast and standardized staining

SAFE BY DESIGN

- Xylene-free, ventilated system
- No dangerous nor toxic reagents used
- Automated to eliminate human error

REPRODUCIBILITY

- Automated, not dependent on operator skills
- Use standard protocols to fit lab requirements
- Pre-program individual lab protocols

HIGH THROUGHPUT

- Handle large cases with multiple specimens
- Simultaneously stain 6 slides at once
- Users can complete other work while staining

FULLY AUTOMATED

- Staining protocols are fully customizable
- Based on pre-set criteria and user protocols
- Ability to create multiple staining protocols

REAGENT MANAGEMENT

- Uses FineFIX for exceptional frozen morphology
- Open system, use your own reagents
- Fully downloadable documentation

UNIQUE FEATURES

8 STEP PROCESSING/STAINING SETUP

PrestoPRO's 10 stations are set up so that you have an 8 step staining protocol with additional vessels for holding, which is ideal for continuous batch loading

MULTIPLE LOADING MODES

Operators can choose either single loading mode or continuous batch loading mode, whichever best fits their workload and case by case requirements

PROCESSING + STAINING

By utilizing Milestone's FineFIX solution, you can now create frozen section slides that mimic formalin fixed, paraffin embedded slides of permanent quality



PrestoPRO utilizes a unique setup for processing and staining, beginning with FineFIX. Preheated at 37°C, FineFIX solution simultaneously performs specimen fixation, dehydration and clearing. Following this step, the eight step staining processing protocols ensure high quality frozen section slides are achieved, every time.

PrestoCHILL

Cryoembedding
of biospecimens



1 MINUTE

Your cryostat


Cutting of sections



~2 MINUTES

PrestoPRO

Automated slide
processing/staining



2-4 MINUTES



US Headquarters

MILESTONE MEDICAL TECHNOLOGIES, INC.

6475 Technology Avenue, Suite F - Kalamazoo, MI 49009

Tel: 269-488-4950 - Toll-free: 866-995-5300 - Fax: 269-488-4949

Email: info@milestonemed.com - www.milestonemed.com