



MILESTONE
H E L P I N G
P A T I E N T S



RoseSTATION

**Mobile FNA Cart,
Arrive Rapidly On-Site**



RoseSTATION

RoseSTATION; the first ever fully equipped and mobile Cytology cart. RoseSTATION can be set up and standardized based on each facilities' standard operating procedures and preferences. When needed, simply unplug it and roll it effortlessly to the desired location. With RoseSTATION, you improve patient care by always arriving rapidly on site, fully prepared and organized for the procedure at hand.

LED LAMP
for additional light

COVERED STORAGE SPACE
for reagents (not included)

SLIDE-OUT DRAWERS
for additional supply storage

UNIVERSAL POWER SUPPLY
for on-board power

BUILT-IN AIR DRYER
for rapid slide drying

SLIDE OUT TRAY
for laptop use

SHOCK ABSORBING PLATFORM
for microscope (not included)

SLIDE OUT DRAWER
for microscope storage



UNIQUE FEATURES

CUSTOMIZABLE SETUP

Each lab is able to set up the RoseSTATION with items required according to their own standard procedures

ON-BOARD POWER SUPPLY

The universal power supply provides up to 4 hours of power. Simply plug it in when not in use

MOBILE

RoseSTATION is always ready to go. Just disconnect it from power and easily push it to the patient's room

ALL-IN-ONE

All procedural items needed are readily available to the operator, assuring rapid turnaround time



US Headquarters

MILESTONE MEDICAL TECHNOLOGIES, INC.

6475 Technology Avenue, Suite F - Kalamazoo, MI 49009

Tel: 269-488-4950 - Toll-free: 866-995-5300 - Fax: 269-488-4949

email: info@milestonemed.com - www.milestonemed.com