




MILESTONE
H E L P I N G
P A T I E N T S



HistoDream EW

Embedding Workstation with
Optional Biospecimen Traceability



HistoDream EW

HistoDream EW is Milestone's alternative to standard embedding workstations. The three independent modules offer flexibility, allowing labs to arrange the embedding workflow in the direction that best suits their daily needs. When paired together with the optional HistoDream TS, tracking system, HistoDream EW allows for visual verification of the cassette contents at grossing compared with the cassette contents at embedding, helping prevent tissue loss.

SET-UP FLEXIBILITY

Set up the thermal and cold modules on either side to best suit your needs



without HistoDream TS

STAND ALONE COLD MODULE

Operates as part of the embedding station or as a stand-alone unit



TRACKING SYSTEM (optional)

View gross images of cassettes with HistoDream TS for full traceability



with HistoDream TS



FULLY EQUIPPED

Programmable heating, paraffin trimmers and heated forceps wells
*Optional heated forceps available

KEY BENEFITS

ADDED FEATURES

Designed with two built-in paraffin trimmers and heated forceps wells for convenience

CASE TRACEABILITY (optional)

Gross image viewing enables detection of mismatched number, size, or specimen type prior to embedding

PROGRAMMABLE HEATING

HistoDream EW is ready at any time, even for urgent or unexpected samples

CAPACITY

Holds up to 300 cassettes, multiple racks, 400 molds, and 5 liters of molten paraffin



US Headquarters

MILESTONE MEDICAL TECHNOLOGIES, INC.

6475 Technology Avenue, Suite F - Kalamazoo, MI 49009

Tel: 269-488-4950 - Toll-free: 866-995-5300 - Fax: 269-488-4949

Email: info@milestonemed.com - www.milestonemed.com