




MILESTONE
H E L P I N G
P A T I E N T S



R-Tracker

Automatic and Efficient
Rack Scanning



R-Tracker

R-Tracker eliminates manual cassette scanning by completely automating this process. With just one-click, R-Tracker simultaneously counts and scans all cassettes and documents each one. R-Tracker detects any missing or unreadable codes, ensuring the integrity of each batch. Complete documentation ensures chain of custody in seconds.

ACCURATE AND RELIABLE

Easily detects cassettes with missing or unreadable barcodes



with built-in PC and monitor

THREE-IN-ONE

Scan, check and track cassettes with double-sensor technology



EASY INTEGRATION

Utilize MileWATCH or your laboratory's LIS system for documentation



without PC and monitor



TIME SAVINGS

Automatically scans all cassettes and their codes in just 30 seconds

KEY BENEFITS

COMPLETE DOCUMENTATION

Automatic scanning provides accurate data on the number of cassettes and corresponding codes

ERROR-FREE SCANNING

A built-in laser scanner easily locates any cassettes with missing or unreadable codes

CASSETTE TRACKING

Integrated with MileWATCH for monitoring and tracking specimens; easily shared through your LIS.

INCREASED PRODUCTIVITY

R-Tracker eliminates manual scanning, a tedious and time consuming task, and frees up tech time



US Headquarters

MILESTONE MEDICAL TECHNOLOGIES, INC.

6475 Technology Avenue, Suite F - Kalamazoo, MI 49009

Tel: 269-488-4950 - Toll-free: 866-995-5300 - Fax: 269-488-4949

email: info@milestonemed.com - www.milestonemed.com