




**MILESTONE**  
H E L P I N G  
P A T I E N T S



# RoseSTATION

Mobile FNA Cart,  
Arrive Rapidly On-Site



# RoseSTATION

RoseSTATION; the first ever fully equipped and mobile Cytology cart. RoseSTATION can be set up and standardized based on each facilities' standard operating procedures and preferences. When needed, simply unplug it and roll it effortlessly to the desired location. With RoseSTATION, you improve patient care by always arriving rapidly on site, fully prepared and organized for the procedure at hand.

**LED LAMP**  
for additional light

**BUILT-IN AIR DRYER**  
for rapid slide drying

**COVERED STORAGE SPACE**  
for reagents (not included)

**SLIDE OUT TRAY**  
for laptop use

**SLIDE-OUT DRAWERS**  
for additional supply storage

**SHOCK ABSORBING PLATFORM**  
for microscope (not included)

**UNIVERSAL POWER SUPPLY**  
for on-board power

**SLIDE OUT DRAWER**  
for microscope storage



## UNIQUE FEATURES

### CUSTOMIZABLE SETUP

Each lab is able to set up the RoseSTATION with items required according to their own standard procedures

### ON-BOARD POWER SUPPLY

The universal power supply provides up to 4 hours of power. Simply plug it in when not in use

### MOBILE

RoseSTATION is always ready to go. Just disconnect it from power and easily push it to the patient's room

### ALL-IN-ONE

All procedural items needed are readily available to the operator, assuring rapid turnaround time



### US Headquarters

MILESTONE MEDICAL TECHNOLOGIES, INC.

6475 Technology Avenue, Suite F - Kalamazoo, MI 49009

Tel: 269-488-4950 - Toll-free: 866-995-5300 - Fax: 269-488-4949

email: [info@milestonemed.com](mailto:info@milestonemed.com) - [www.milestonemed.com](http://www.milestonemed.com)