




MILESTONE
H E L P I N G
P A T I E N T S



UltraGROSS



The New Benchmark
in Grossing Stations



Customizable Grossing Station

UltraGROSS

UltraGROSS addresses the needs of the lab of the future with its focus on modularity, digital documentation, and operator safety. With UltraGROSS, the Pathologist's Assistant can generate reports that provide complete documentation and traceability from the time the specimen is grossed. Ergonomically designed with enhanced user safety features, UltraGROSS also offers fully customizable user configurations to fit the needs of any lab, today and tomorrow.






SEE ULTRAGROSS IN ACTION

** The unit pictured is the UltraGROSS Virchow, which is the only available UltraGROSS model capable of fitting 2 users, as shown*

AVAILABLE MODELS

Choose one of our three pre-set configurations or request a tailor-made unit to fit your specific requirements. Whatever your specifications, there's an UltraGROSS unit to fit your needs!

	GOLGI <i>Edition</i>	ROSAI <i>Edition</i>	VIRCHOW <i>Edition</i>
			
Available Workstation Length	150cm or 180cm	150cm or 180cm	180cm
Included Digital Documentation	MacroPATH camera x1 (no CPU)	MacroPATH System x1 (with CPU)	MacroPATH System x2 (with CPU)
Number of Operators	1	1	2

KEY BENEFITS

ENHANCED USER SAFETY

- Horizontal sliding safety shield
- Powerful downdraft and backdraft exhaust
- Designed to eliminate formalin exposure

TAILOR-MADE MODULARITY

- Design around current and future needs
- Moveable modules across the entire basin
- Interchangeable layouts to fit user needs

SMART USER INTERFACE

- 8" smart touchscreen terminal
- Capability to upload external SOP files
- All user functions at your fingertips

DIGITAL DOCUMENTATION

- Macro digital documentation with MacroPATH
- Ensures complete traceability/chain of custody
- Enhances quality of diagnostic reporting

OPTIMAL ERGONOMICS

- Easily height adjustable working area
- Casters with safety locks for mobility
- Set-up flexibility to fit specific operator needs

FUTURE FLEXIBILITY


- Easily modified if future needs change
- Moveable sink basin. Right, left, or center
- Full set-up flexibility, no need for technician

UNIQUE FEATURES



MACRO DIGITAL DOCUMENTATION

Full documentation is assured via the built-in MacroPATH macro digital imaging system (optional), the most widely installed imaging system in the world. MacroPATH enhances the quality of diagnostic reports for tumor boards, and allows for specimen tracking downstream.



OPTIMAL ERGONOMICS

The working area height can be easily adapted to the position most comfortable for each individual user through the on-board, user-friendly, touchscreen control terminal.



ENHANCED SAFETY

Built for safety with spent and fresh formalin tanks with safety valves, built-in backdraft and downdraft ventilation, and a sensor to detect operator presence at the bench.



SMART USER INTERFACE

An 8" industrial touchscreen terminal allows the operator to have full control of the unit, including the exhaust fan, the height of the unit, the light, and the operator sensor.

RELATED PRODUCTS



MacroPATH:
Digital Imaging System



MileWATCH:
24/7 Monitoring & Tracking System



T-Tracker:
Automatic Rack Scanner



HistoDream TS:
Biospecimen Tracking System



US Headquarters

MILESTONE MEDICAL TECHNOLOGIES, INC.

6475 Technology Avenue, Suite F - Kalamazoo, MI 49009

Tel: 269-488-4950 - Toll-free: 866-995-5300 - Fax: 269-488-4949

email: info@milestonemed.com - www.milestonemed.com