




MILESTONE
H E L P I N G
P A T I E N T S



V-Tracker


Digital Vial Imaging System
for Specimen Documentation



Automatic Vial Scanner


V-Tracker

Milestone's V-Tracker is an automatic imaging tool designed to simplify and improve specimen accessioning. V-Tracker's software captures information from the specimen container, including barcodes, QR codes, and any additional notes present on the container. It then consolidates all this information into a single image, ensuring comprehensive documentation of each vial scanned. Taking it one step further, V-Tracker can capture an image of the inside of the specimen container for more advanced traceability. Utilizing standardized imaging with V-Tracker, labs can enhance accuracy and reliability, and allow for seamless specimen tracking down stream.



CAPTURE EVERY DETAIL

VIAL IN AND
SCANNED



INFORMATION
ACQUIRED

360° scanning enables the capture of the entire specimen container, including:

- The container label
- All codes present on the container
- Any notes or additional information

KEY BENEFITS

OPTIMIZED WORKFLOW

Streamlines the specimen acceptance and tracking process for error-free operation and efficiency

COMPLETE DOCUMENTATION

Automatically prioritizes the master code and records all other codes as secondary for traceability

STANDARDIZED TRACKING

Scan all specimen vials in the same way, every time, to ensure data collection standardization

HIGH FLEXIBILITY

Scans specimen containers of all different shapes and sizes, as well as large specimen containers



US Headquarters

MILESTONE MEDICAL TECHNOLOGIES, INC.

6475 Technology Avenue, Suite F - Kalamazoo, MI 49009

Tel: 269-488-4950 - Toll-free: 866-995-5300 - Fax: 269-488-4949

Email: info@milestonemed.com - www.milestonemed.com